


ISOTOPE-LABELED PROTEINS FOR NMR

- ▶ High-yield production of isotope-labeled protein for NMR
- ▶ Optimised expression protocols for prokaryotic and eukaryotic expression
- ▶ Bespoke labeling strategies to facilitate specific NMR assays
- ▶ State of the art purification and protein analysis
- ▶ Short turn-around times to obtain high-quality, labeled protein at low cost
- ▶ Over ten years experience in isotope-labeling & a track record of >100 labeled proteins



IN-HOUSE OPTIMISED CULTURE MEDIA AND PROTOCOLS FOR COST-EFFICIENT EXPRESSION OF ISOTOPE LABELED PROTEINS:

	<i>E. coli</i>	<i>P.pastoris</i>	Insect cells	Mammalian cells
Ubiquitous labeling	¹⁵ N, ¹³ C, ² H	¹⁵ N, ¹³ C	¹⁵ N, ¹³ C-labeling possible, but affiliated with high media costs	
Site-specific labeling*	¹³ C- Val, Leu/ Ile using α-keto -isovalerate/ -butyrate	–	¹⁵ N- and/or ¹³ C for: Val, Leu, Met, Tyr, Phe, Ile, His, Trp, Thr	–
Media properties	NH ₄ as N-source, glucose (<i>E. coli</i>) and methanol (<i>P.pastoris</i>) as C-source, dedicated trace metal and vitamin mixes		taylor-made minimal media	–

Evotec Protein Production Services

▶ *Highly experienced team and state-of-the-art technology platform for membrane protein production, isotope-labeled protein preparation, high throughput protein purification and structural biology services*


▶ *Direct communication with project scientists and full access to raw data*

▶ *Reagent production services can be accessed as a standalone function or as integral part of wider collaborative programmes*

▶ *Links to key academic groups from Oxford University & Diamond light source, UK*

CASE STUDY 1

- ▶ Expression of a ¹⁵N-labeled serine protease in *Pichia pastoris* was optimised, e.g. expression time/temperature
- ▶ The purified protease was submitted to an NMR fragment screen using ¹⁵N-HSQC
- ▶ Active-site ¹H-¹⁵N-resonances could be identified by analysing chemical shift perturbations induced by benzamidine binding:



CASE STUDY 2

- ▶ Expression of a ¹⁵N/¹³C-Met/Leu-labeled kinase in insect cells
- ▶ Due to the presence of a unique Leu-Met dipeptide, the ¹H-¹⁵N-signal of the gatekeeper Met could be identified using CCLP-HSQC (see below). This facilitated compound binding studies where X-ray crystallography failed

